



## Technical specification

**Model:** DW- 2 5YL270

**Effective Capacity:** 270 Liters

**Temperature range:** -10°C ~ -25°C

### Feature

- Vertical cabinet (upright)
- Temperature range: from -10°C.to -25°C, can be set freely.
- microprocessor temperature controller
- LCD display, to make clear the temperature information
- ABS material interior, Outer body is painted steel board.
- the safety lock on the door
- Rotatable sure-grip door handle
- Perfect audible/visual alarm: High or low temperature alarm, system failure alarm
- 90 mm extra thick insulation
- The Danfoss compressor (1unit) and EBM fan motor.
- Password protected configuration page
- 4 Casters provided for easy handling
- 7units drawers made of plastic inside.
- Refrigerant as R507, CFC free
- Power supply: 220V /50Hz 1 phase, can change as 220V 60HZ or 110V 50HZ or 110V 60HZ
- Certificate: ISO9001, ISO14001, ISO13485
- Optional: chart recorder

### Specifications

Model	capacity	External size (W*D*H)mm	Inside size(like (W*D*H)mm	Input power	Weight (Nt/Gt)
DW-25YL270	270liters	623*544*1680	438*388*1580	370W	65/70