

Model: DW-YW110A/166A/196A/226A/358A/508A, Chest Type Freezer**Effective capacity: 110/166/196/226/358/505 Liters. Temperature range: -10°C ~ -25°C****Features:**

- Chest type.
- Inside temperature from -10°C to -25°C, can be set freely.
- Digital temperature display.
- Round corner design of doors, with lock.
- Outer body is painted steel board, inside is aluminum.
- 40mm foaming insulation.
- China good quality compressor (1 unit).
- 1 unit basket made of steel wire are convenient to store articles.
- Delayed start and safe stop interval between restart and being terminated.
- Quick freezing function, 2-3 hours come to -25°C very fast.
- 4 units Casters for easy handling.
- Refrigerant as R134a, CFC free.
- Certificate: ISO9001, ISO14001, ISO13485.
- Power supply: 220V/150Hz 1 phase, can change as 220V/60Hz or 110V/50Hz or 110V/60Hz.
- **Optional: Basket.**

Model	Capacity	External size(mm)	Inside size(mm)	Input power	Weight (N.T/G.T)
DW-25YW110A	110 Liters	W549xD549xH845	W410xD410xH654	145W	31Kg/46Kg
DW-25YW166A	166 Liters	W556xD906xH937	W430xD780xH480	160W	46Kg/71Kg
DW-25YW196A	196 Liters	W556xD1056xH937	W430xD930xH480	180W	51Kg/67Kg
DW-25YW226A	226 Liters	W556xD1206xH937	W430xD1080xH480	207W	56Kg/71Kg
DW-25YW358A	358 Liters	W730xD1204xH968	W530xD1080xH625	320W	82Kg/97Kg
DW-25YW508A	508 Liters	W730xD1554xH968	W530xD1400xH685	375W	100Kg/115Kg

**Model: DW-40FL90/135/270, DW-25YL270, Upright Type Freezer****Effective Capacity: 90/135/270 Liter.****Temperature range: -10°C ~ -40°C, -10°C ~ -25°C****Features:**

- Vertical cabinet (upright).
- Temperature range: from -10°C to -40°C, -10°C to -25°C, can be set freely.
- Microprocessor temperature controller.
- LCD display, to make clear the temperature information.
- Elegant orange plastic coating.
- ABS interior.
- The safety lock on the door.
- Rotatable sure-grip door handle.
- Perfect audible/visual alarm: High or low temp. alarm, system failure alarm.
- 90 mm extra thick insulation.
- The Danfoss compressor (1unit) and EBM fan motor.
- Password protected configuration page.
- 4 Casters provided for easy handling.
- 4units drawers(plastic) for DW-FL90, 5units drawers for DW-FL135, 7units drawers for DW-FL270.
- Refrigerant as R507, CFC free.
- Power supply: 220V/50Hz 1 phase, can change as 220V/60Hz or 110V/50Hz or 110V/60Hz.
- Certificate: CE mark for DW-FL90/270, ISO9001, ISO14001, ISO13485.
- **Optional: chart recorder, Italy Carel controller.**

Model	Capacity (Liter)	External size(mm)	Inside size(mm)	Input power	Weight (N.T/G.T)
DW-40FL90	90	W623xD544xH842	W438xD388xH539	280W	44Kg/48Kg
DW-40FL135	135	W623xD544xH1123	W438xD388xH794	280W	54Kg/60Kg
DW-40FL270	270	W623xD544xH1680	W438xD388xH1580	370W	65Kg/70Kg
DW-25YL270	270	W623xD544xH1680	W438xD388xH1580	370W	65Kg/70Kg

-86°C Models	Rack Size (mm)	Racks/Freezers	Boxes/Racks	Total Boxes / Freezer
DW-86HL218	140 x 435 x 336	6	20	120
DW-86HL328	140 x 435 x 280	12	20	240
DW-86HL388	140 x 550 x 225	12	20	240
DW-86HL538	140 x 550 x 225	16	20	320
DW-86HL688	133 x 510 x 133	20	20	400
DW-86HL768	133 x 510 x 133	24	20	480
DW-86HL828	140 x 680 x 320	24	30	720

-65°C Models	Rack Size (mm)	Racks/Freezers	Boxes/Racks	Total Boxes / Freezer
DW-65GL218	140 x 435 x 336	6	20	120
DW-65GL328	140 x 435 x 280	12	20	240
DW-65GL388	140 x 550 x 225	12	20	240
DW-65GL538	140 x 550 x 225	16	20	320
DW-65GL828	140 x 680 x 320	24	30	720

-86°C Models	Rack Size (mm)	Racks/Freezers	Boxes/Racks	Total Boxes / Freezer
DW-86HW138	143 x 140 x 540	8	16	128
DW-86HW328	143 x 140 x 560	24	16	384

-65°C Models	Rack Size (mm)	Racks/Freezers	Boxes/Racks	Total Boxes / Freezer
DW-65GW138	143 x 140 x 540	8	16	128
DW-65GW328	143 x 140 x 560	24	16	384

-40°C Models	Rack Size (mm)	Racks/Freezers	Boxes/Racks	Total Boxes / Freezer
DW-40FL531	140 x 550 x 225	16	20	320

Option: SMS Temperature & Humidity Alarm

AL-SMS-THR Series, Smart Energy Control ... Energy Saving, Environment & Properties Protection



3 simple steps starting the unit:

- Insert the SIM card
- Plug in the sensors connectors
- Connect the AC power cord.

Specifications:

- AC 90~260V Auto Select
- Internal 2Ah Rechargeable Battery
- Non-Stop Operation in AC power loss
- SMS Alert on AC power loss.

Application:

Temperature & Humidity SMS Alert Controller can be utilized in monitoring the environmental condition for commercial, medical, food , storage, weather station, military, manufacturing industries.

It makes use of SMS to alert users of temperature and humidity conditions via SMS. By implementing the bundled Server/Client Software and GSM modem, a management system of multi-controllers can be built up to monitor a number of remote sites in one single centre.

Features:

- Alarm Alert Modes—SMS, Phone Dial or Both.
- Each alarm is associated with its own alarm SMS text.
- Alarm SMS text is programmable by user remotely via SMS or Setup Software.
- Independent text message for each alarm . Independent text message for close or open triggered alarm.
- Configuration can be modified by user remotely via SMS or Setup Software.
- Arm or Disarm by SMS or Setup Software.
- Relay triggered by alarm or SMS command.
- Verification of command by return SMS.