



DENT-P15E, Portable Dental Equipment Dental Unit

The portable dental unit is dedicated to provide convenient and professionally stable services for dentists to go out for medical treatment with portable dental equipment. Professional and high-quality medical portable dental equipment wins praise from doctors and patients. It can be used as a backup machine for dental clinics or as a special machine for veterinary treatment.

Portable dental unit is a unit of dental chair basic functions that we provide for outpatient dental treatment and meets the needs of the clinic. Convenient, durable, professional, etc., which is an indispensable partner for dental professionals

Luggage case portable dental unit with Lightest Weight (**15KG only**)
 Compatible with 110V & 220V voltage
 Available with **portable power source**.

ABS plastic has high strength, light weight, high surface hardness, very smooth, easy to clean, stable size, good creep resistance, good electrical insulation, almost no influence on temperature, Humidity and frequency. It is used in most environments. Molded ABS material suitcase design make the self-contained portable dental unit more light weight, pressure resistant, easy to transport.

The portable dental equipment control by a foot pedal which can be arranged inside the self-contained portable dental unit. Multifunctional two-in one foot pedal which control both of the hand piece and scaler, more convenient.

Applications:

Dental Hospital, Dental Clinics, Army and Outdoor Dental Care or Treatment Establishment, etc.

Portable dental unit is the most convenient portable dental equipment for mobile dental clinic. The function of dental chair is concentrated in a small box, which saves a lot of space for mobile dental van or mobile dental office, and it is easy to pull the portable dental unit for mobile dental services. Especially for dentist home visit. And because of portable dental unit be with portable dental compressor & dental handpiece, dental scaler, it can be used as portable dental delivery unit or portable dental suction unit.

- Lightweight easy to carry and beautiful appearance design makes it more convenient for outpatient treatment. The professional treatment equipment makes the patient feel more at ease.
- A more concise external construction system reduces the possibility of breakdown and enables the equipment to be stable and high quality for a long time
- Taking into account the actual work needs, improving the innovative and humanized details of design, providing more intimate services & comfortable working procedures for medical personnel
- It has high volume external circulation suction incorporating clean airflow, maintaining optimum suction power when used in conjunction with the hand piece.
- Compact structure consistent with doctor's habits to make them operate more easier
- Using when install the parts correctly and connect with power. And the unit controlled with foot control pad.
- Oil free motor without any lubricated oil, drainage bottle is more convenient to remove and disinfect. **Good quality lies in the details**



Model	DENT-P15E
Power supply	100-240V, 160W
Weight	15KG
Max. pressure	0.8 Mpa
Noise	35dB