

# Fraction Collector



## FRC-100/160, Fraction Collector

### Features:

- Keebio fluid stop valve, test-tube disc resistant to organic solvents (optional)
- 35°bevel control panel design. There is no need to bend down to easily control.
- Super blue film LCD super clear display. new microcomputer chip control.
- Unique collection plate digital marking design. each tube corresponds to a number, to ensure that the tube is accurately marked.
- Collection function: automatic collection, arbitrary selection, no deviation and accurate positioning when tube changing.
- Timing function: collect time, invert, clockwise timing mode is optional. Free conversion of tube with same specification.

Model	FRC-100	FRC-160
<b>Tube quantity</b>	100	160
<b>Volume</b>	12ml	5ml
<b>Regular collection range</b>	1s - 199 hours, 59 minutes & 59 seconds, optional	
<b>Drop collection range</b>	1 drop - 9999 drops, optional	
<b>Peak collection range</b>	0-100mv	
<b>Display mode</b>	LCD display	
<b>Tube Plate</b>	The tube plate is marked with a circular number	
<b>Power off data storage</b>	Decade	
<b>Continuous working time</b>	More than 200 hours	