


ROL-01U, Mini Roller



ROL-01U



Stacking Bars for
Stackable Mini
Roller, 4 per pack

Model	ROL-01U
Number of Roller Bars	6 Pcs Standard, 5 Positions
Speed	0.5~80rpm
Roller Dimensions (DxL)	φ28x269mm
Available Tubes or Bottles	Max. φ120mm, bottles/tubes
Max. Load	6.5kg, with container, 3x1500ml bottles
Digital Display	Roller Speed Display
Timing Range	0~999min or continuous operation
Included Accessory Standard in Packaging	MRA-02, MRA-03, MRA-05
Ambient Temperature	5~60°C
Relative humidity	≤95%
Max. Power	10W
Dimensions (WxDxH)	378x360x72mm
Net Weight	6.5kg
Stacking	Optional