



Laboratory Equipment Manufacturer
www.mrclab.com



Features

- 1.100w brushless dc motor,single chip microcomputer control
- 2.High-definition 4.3" display
- 3.Supporting wafer from 1cm-150mm(6")
- 4.Pre-set single or multi-step coating operate
- 5.Prevent glue blocking design
- 6.PP anti- corrosion chamber, replaceable
- 7.HDPE vacuum adhesive disc

Model No.	SPC-60
Display	4.3" LCD Touch Panel
Sample range	10-150mm
Speed resolution	±1RPM
Rotating speed	100~12000RPM

Accelerated speed	100~3000rpm/s
Single step timing	0.1~3000s
Control	10 sets programmable control,10steps
Motor power	Brushless DC Motor 100w
Vacuum pump	Oiless vacuum pump 15L/min
Working voltage	AC 220V, 50Hz,5A
Dimension	358*256*258mm



3, Hagavish st. Israel 58817 Tel: 972 3 5595252, Fax: 972 3 5594529

mrc@mrclab.com

MRC. 5.19