



TRO-1000N, Tube Roller, 16 Rollers

It can accommodate up to 16×5~10ml blood vessels, or 6×50ml tubes, with different silicon mattress accessories.

Model	TRO-1000N
Tilt angle	30±3°
Load capacity	1kg
Motor type	Brushless DC motor
Speed range (rpm)	0~30
Operation mode	Continuous
Dimension (W×D ×H)	405×160×120mm
Voltage	100~240V, 50/60Hz
Power	15W
Weight	2kg
Permissible ambient temperature & humidity	5~40°C, 80%