




197000-11, PES waste bottle

| | | | |
|--|-------------------------|--------------------------|--------|
|  <p>197000-11</p> | Material Specifications | Capacity | 1200ml |
| | | Port thread(I.D. tubing) | 8mm |
| | | Graduated increment | 100ml |
| | | Spill-proof buoy | Yes |
| | | Waste bottle cover | Yes |
| | | Compatible Stopper | No.16 |
| Material | Container | PES | |
| | Cover | PES | |
| | Spill-proof buoy | PP | |
| | Stopper | Silicon | |


197000-12, PES waste bottle

| | | | |
|--|-------------------------|--------------------------|--------|
|  <p>197000-12</p> | Material Specifications | Capacity | 1200ml |
| | | Port thread(I.D. tubing) | 5mm |
| | | Graduated increment | 100ml |
| | | Spill-proof buoy | Yes |
| | | Container | PES |
| | | Cover | PES |
| Material | Spill-proof buoy | PP | |


197000-13, Suction / Pressure bottle

| | | | |
|--|-------------------------|--------------------------|--------|
|  <p>197000-13</p> | Material Specifications | Capacity | 1200ml |
| | | Port thread(I.D. tubing) | 5mm |
| | | Graduated increment | 100ml |
| | | Overpressure protection | Yes |
| | | Container | PES |
| Material | Cover | PES | |


167200-33, PC waste bottle

| | | | |
|--|-------------------------|--------------------------|--------|
|  <p>167200-33</p> | Material Specifications | Capacity | 3000ml |
| | | Port thread(I.D. tubing) | 8mm |
| | | Graduated increment | 100ml |
| | | Spill-proof buoy | Yes |
| | | Container | PC |
| | | Cover | PC |
| Material | Spill-proof buoy | PP | |


197000-36, PP funnel lid

| | | | |
|--|-------------------------|-------------------|------------------|
|  <p>197000-36</p> | Material Specifications | Compatible funnel | 300ml PES funnel |
| | | Vent | 3 |
| | | Vent Fitting | Female luer slip |
| | | Syringe filter | Ø25mm, 0.2µm |
| | | Body | PP |
| | | Cap | Silicon |
| | | Syringe filter | PTFE |


167200-34, Glass waste bottle

| | | | |
|--|-------------------------|--------------------------|--------|
|  <p>167200-34</p> | Material Specifications | Capacity | 4000ml |
| | | Port thread(I.D. tubing) | 8mm |
| | | Graduated increment | 100ml |
| | | Spill-proof buoy | Yes |
| | | Container | Glass |
| | | Cover | S.S |
| Material | Spill-proof buoy | PP | |


167100-20/52, Membrane filters

| | | | |
|---|---------------------------------|-----------------|--------|
|  <p>197100-20</p> | 167100-20 Glass fiber filter | Brand | Pall |
| | | Filter diameter | 47 mm |
| | | Pore size | 1µm |
| | 167100-52 MCE filter | Brand | Pall |
| | | Filter diameter | 47 mm |
| | | Pore size | 0.45µm |
| Material | Material | MCE | |


167100-13, Glass flask with drain

| | | | |
|--|-------------------------|--------------------------|--------------------|
|  <p>167100-13</p> | Material Specifications | Capacity | 1000ml |
| | | Port thread(I.D. tubing) | 8mm |
| | | Graduated increment | 100ml |
| | | Spill-proof buoy | Yes |
| | | Compatible Stopper | No.16 |
| | | Container | Borosilicate glass |
| Material | Connector | PP | |
| | Stopper | Silicon | |

194225-01,167200-35/36, Syringe filters

| | | | |
|---|--------------------------------|------------------------|------------|
|  <p>167200-35</p> | 194225-01 Syringe filter | Brand | Pall |
| | | Filter dia./ Pore size | 25mm/0.2µm |
| | | Shell diameter | 30mm |
| | 167200-35 Small disc filter | Filter material | PTFE |
| | | Shell material | PP |
| | | Brand | Pall |
| | 167200-36 Large disc filter | Filter dia./ Pore size | 35mm/0.2µm |
| | | Shell diameter | 45mm |
| | | Filter material | PTFE |
| Material | Shell material | PP | |
| | Brand | Pall | |
| | Filter dia./ Pore size | 50mm/0.2µm | |
| Material | Shell diameter | 72mm | |
| | Filter material | PTFE | |
| Material | Shell material | PP | |

167200-13, Glass flask with drain

| | | | |
|--|-------------------------|--------------------------|--------------------|
|  <p>167200-13</p> | Material Specifications | Capacity | 2000ml |
| | | Port thread(I.D. tubing) | 8mm |
| | | Graduated increment | 200ml |
| | | Spill-proof buoy | Yes |
| | | Compatible Stopper | No.18 |
| | | Container | Borosilicate glass |
| Material | Connector | PP | |
| | Stopper | Silicon | |

PES: Polyethersulfone, PC: Polycarbonate, PP: Polypropylene