




197000-11, PES waste bottle

 <p>197000-11</p>	Material Specifications	Capacity	1200ml
		Port thread(I.D. tubing)	8mm
		Graduated increment	100ml
		Spill-proof buoy	Yes
		Waste bottle cover	Yes
		Compatible Stopper	No.16
Material	Container	PES	
	Cover	PES	
	Spill-proof buoy	PP	
	Stopper	Silicon	


197000-12, PES waste bottle

 <p>197000-12</p>	Material Specifications	Capacity	1200ml
		Port thread(I.D. tubing)	5mm
		Graduated increment	100ml
		Spill-proof buoy	Yes
		Container	PES
		Cover	PES
Material	Spill-proof buoy	PP	


197000-13, Suction / Pressure bottle

 <p>197000-13</p>	Material Specifications	Capacity	1200ml
		Port thread(I.D. tubing)	5mm
		Graduated increment	100ml
		Overpressure protection	Yes
		Container	PES
Material	Cover	PES	


167200-33, PC waste bottle

 <p>167200-33</p>	Material Specifications	Capacity	3000ml
		Port thread(I.D. tubing)	8mm
		Graduated increment	100ml
		Spill-proof buoy	Yes
		Container	PC
		Cover	PC
Material	Spill-proof buoy	PP	


197000-36, PP funnel lid

 <p>197000-36</p>	Material Specifications	Compatible funnel	300ml PES funnel
		Vent	3
		Vent Fitting	Female luer slip
		Syringe filter	Ø25mm, 0.2µm
		Body	PP
		Cap	Silicon
		Syringe filter	PTFE


167200-34, Glass waste bottle

 <p>167200-34</p>	Material Specifications	Capacity	4000ml
		Port thread(I.D. tubing)	8mm
		Graduated increment	100ml
		Spill-proof buoy	Yes
		Container	Glass
		Cover	S.S
Material	Spill-proof buoy	PP	


167100-20/52, Membrane filters

 <p>197100-20</p>	167100-20 Glass fiber filter	Brand	Pall
		Filter diameter	47 mm
		Pore size	1µm
	167100-52 MCE filter	Brand	Pall
		Filter diameter	47 mm
		Pore size	0.45µm
Material	Material	MCE	


167100-13, Glass flask with drain

 <p>167100-13</p>	Material Specifications	Capacity	1000ml
		Port thread(I.D. tubing)	8mm
		Graduated increment	100ml
		Spill-proof buoy	Yes
		Compatible Stopper	No.16
		Container	Borosilicate glass
Material	Connector	PP	
	Stopper	Silicon	

194225-01,167200-35/36, Syringe filters

 <p>167200-35</p>	194225-01 Syringe filter	Brand	Pall
		Filter dia./ Pore size	25mm/0.2µm
		Shell diameter	30mm
	167200-35 Small disc filter	Filter material	PTFE
		Shell material	PP
		Brand	Pall
	167200-36 Large disc filter	Filter dia./ Pore size	35mm/0.2µm
		Shell diameter	45mm
		Filter material	PTFE
		Shell material	PP
		Brand	Pall
		Filter dia./ Pore size	50mm/0.2µm
Material	Shell diameter	72mm	
	Filter material	PTFE	
Material	Shell material	PP	

167200-13, Glass flask with drain

 <p>167200-13</p>	Material Specifications	Capacity	2000ml
		Port thread(I.D. tubing)	8mm
		Graduated increment	200ml
		Spill-proof buoy	Yes
		Compatible Stopper	No.18
		Container	Borosilicate glass
Material	Connector	PP	
	Stopper	Silicon	

PES: Polyethersulfone, PC: Polycarbonate, PP: Polypropylene