

GL45 Versatile Cap

Applications



Waste Bottle

Suitable for collecting waste media by using with GL 45 bottle in connection with suction pump and accessories.




Moisture Trap Bottle

Available to absorb moisture with desiccant which turns pink when saturated.




Vapor Absorption (Neutralization) Bottle

Ideal for absorbing acid or alkaline vapor to prevent harmful vapors from spreading into the environment or pump.



Transfer bottle

Ideal for transferring sterile media stored in GL 45 bottle to another container or reactor with peristaltic pump in a closed system.



Safety Bottle

Works as a safety bottle to prevent liquid from being sucked into pump.



Ordering Information

197000-60	Suction Cap for Storage Bottle (1/pk)
197000-60x10	Suction Cap for Storage Bottle (10/pk)
199100-68	Silicone Tube (Φ 4.2x Φ 8 mm), 200cm
197000-60-500	500 ml Storage Bottle with suction cap
197000-60-1000	1000 ml Storage Bottle with suction cap