



BOD-80E, BIOCHEMISTRY INCUBATOR 80 LITER

Features:

- ▶ LED display Controller.
- ▶ LCD display.
- ▶ Polished stainless steel chamber.
- ▶ Round angle structure , easy to clean.
- ▶ Universal castors for easy movement.
- ▶ CFC-free eco-friendly refrigerant.
- ▶ Universal Castors for easy movement.
- ▶ CFC-free environmental friendly refrigerant.
- ▶ Microprocessor PID temperature control, accurate and reliable.
- ▶ Power socket integrated inside for additional small equipment.

Safety Function:

- ▶ Temperature protection threshold value can be set.
- ▶ Auto-off protection and alarm for abnormal situation.
- ▶ Parameters can be saved automatically when power off or system error.

Optional Accessories:

- ▶ Test port.



INSIDE

Technical Parameters:

MODEL	BOD-80E
Capacity	80 Liter
Temperature Range	0°C ~ 65°C
Temperature Fluctuation	±0.5°C
Temperature Uniformity	±2°C
Time Range	0-9999min
Refrigerant	R134a (CFC free, eco-friendly)
Standard Accessory	BOD socket
Consumption	180W
Power Supply	AC220V±10%, 50/60Hz
Internal Size(mm)	400x370x550
External Size(mm)	540x570x1070
Package Size(mm)	700x730x1260
Gross Weight(kg)	85