

### BWLC-TB Series, Precision Balances





BWLC-TB balances are designed for fast and precise determination of mass in laboratory conditions.

They can be used in locations where no power (230V) is accessible as they are equipped with internal rechargeable battery.

Balances are equipped with stainless steel weighing pan, RS 232 output and backlit LCD display.

#### Optional accessories:

- Printers
- Table for scale
- Computer software
- RS 232/RS 485 Converter
- External power output



Model	BWLC-200-TB	BWLC-2000-TB
Max capacity	200g	2000g
Readability	0.001g	0.01g
Tare range	-200g	-2000g
Repeatability	0.002g	0.01g
Linearity	±0.004g	±0.03g
Pan size	Ø100 mm	128 x 128 mm
Working temperature	+15° - +30°C	
Power supply	230V / 11V AC or 120V / 11V AC, 8 6 AA NiMH accumulator	
IP rating	IP 43	
Display	LCD (backlit)	
Net / Gross weight	1.3/2kg	
Packaging size	330x220x140mm	