





Scales WPT-6-15DC and WPT 10/20DC are precise electronic weighing equipment manufactured in class III and designed for weighing infants. The weighing pan is made of Owispol polystyrene and the display is in backlit LCD version.

Each scale is powered by mains (230V/11V AC). The indicator has also container for batteries (6xAA - an alternative Power source for the scale).

Measuring indicator Type-31 is standard equipment of this scale.

**Optional accessories:**

- Thermal printer
- Computer software



Model	WPT-6-15D	WPT-10-20D
Maximal capacity	6/15kg	10/20kg
Readability	2/5g	5/10g
Tare range	-15kg	-20kg
Working temperature	0 ÷ +40°C	
Interface	RS 232	
Power supply	230V/11V AC or 120V/11V AC and 6xAA NiMH accumulator	
Display	LCD (backlit)	