

Scales designed for weighing patients in clinics. Using height-meter, it is possible to measure the height of a patient up to 2 meters. Scales are equipped with backlit LCD display, Power supply either from mains or from batteries.



Personal scales are manufactured in two versions:

- with height-meter (indicator installed on the pillar of height - meter).
- without height-meter (indicator attached to the weighing platform).

Measuring indicator Type-31 is standard equipment of this scale.

**Optional accessories:**

- Thermal printer
- Computer software



Model	WPT-60-150OW WPT-60-150O	WPT-100/200OW WPT-100/200O
Maximal capacity	60/150kg	100/200kg
Readability	20/50g	50/100g
Tare range	-150kg	-200kg
Working temperature	0 ÷ +40°C	
Output signal	RS 232	
Power supply	230V / 11V AC or 120V / 11V AC and 6×AA NiMH accumulator	
Display	LCD (backlit)	